




SIRIS SPIRAL SPLINT (SP)

Thumb IP Hyperextension or Trigger Thumb

- 1 Choose an EZ-Sizer that, when rotated on its side, is snug when it slides over the IP joint.




- 2 Measure the length from mid-phalanx to mid-phalanx.



- 3 *The spiral splint can be worn to block hyperextension and provide mild medial/lateral stability to the joint.*



The spiral splint can be worn for trigger thumb by positioning the center band to cross diagonally over the dorsal surface of the IP joint.



SPLINT ADJUSTMENTS AND FITTING:

- Splints can be adjusted for circumferential fit by either tightening or opening the spiral.
- Splint length can be adjusted by spreading the ends apart or squeezing the ends towards each other.

SIRIS™ MEASURING INSTRUCTIONS

SILVER RING SPLINT COMPANY | PHONE: 434.971.4052 or 800.311.7028 | FAX: 434.971.8828 or 888.456.8828

CustomerService@SilverRingSplint.com | www.SilverRingSplint.com