



Note: The splint size on the inside of each ring does not indicate the actual size of the splint. Please refer to sizes listed on invoice.

1. The proximal ring, the large end, goes on the finger first.



2. Insert the finger into the proximal ring. Flex the finger to align the top of the finger with the distal ring.



3. The finger must be fully flexed while sliding the proximal ring over the PIP joint. It might be tight but once over the joint it should be comfortable.



4. Gently push the splint back along the finger until the splint is in place.



5. The angle between the rings may be adjusted to get the desired fit and degrees of flexion.

