



# SIRIS EDS SPLINT (EDS)

*Hyperextension of MP Joints*

---


## SPECIAL ORDER SPLINT

---


- 1** Choose a Round Ringsizer that slides easily over the PIP joint and can be positioned at the base of the finger.



- 2** With the Ringsizer still in place, measure the length from the distal palmar crease to the distal edge of the Ringsizer.



- 3** With the Ringsizer still in place and with the finger in slight flexion, measure the length from the center of the MCP joint to the distal edge of the Ringsizer and add 2cm to that measurement.



---

### SIRIS™ MEASURING INSTRUCTIONS

SILVER RING SPLINT COMPANY | PHONE: 434.971.4052 or 800.311.7028 | FAX: 434.971.8828 or 888.456.8828

*CustomerService@SilverRingSplint.com | www.SilverRingSplint.com*