


SIRIS™

NEW THUMB SOLUTIONS!

SIRIS™ STABLE-THUMB

Stabilizes the MP joint in extension (joint must be passively correctable)

- For Z Deformity of the thumb (Type I and Type II thumb deformity)
- Extends the thumb MP joint which reduces IP joint hyperextension
- Increases power and prehension



Low profile, cosmetic splints are essential for people with chronic hand problems allowing them to live gracefully.

SIRIS™ STABLE-THUMB WITH PVX (PROXIMAL VOLAR EXTENSION)

Prevents MP motion and reduces hyperextension

- For swan neck deformity and adduction at MP joint (Type III thumb deformity)
- Frequently used for OA and RA when there are CMC joint changes (pre or post CMC arthroplasty with an adducted prehension pattern)
- Proximal volar extension supports under the head of the metacarpal joint to stabilize the MP joint and reduce hyperextension
- Increases pinch force, improves prehension and decreases pain



- To measure, find the smallest round ring size which will go over the IP joint
- Specify PVX (Proximal Volar Extension) if needed
- Bracelet comes standard; no measurement needed for PVX length

P.O. Box 2856
Charlottesville, VA 22902
Fax: 888-456-8828
Web: www.silverringsplint.com

Silver Ring Splint Company

800-311-7028