



SIRIS 2¹/₂ BUDDY RING SPLINT (2¹/₂ BR)


Adduction of Fingers, Prevent Scissoring or Aligns and Supports MP Joints

- 1** Choose Round Ringsizers that fit at the base of the three adjacent fingers being splinted: index, middle and ring fingers or middle, ring and little fingers.



(Actual size shown)

- 2** Specify band width: 3 or 4mm.



- 3** *Note: The half ring can be worn under the finger to correct an extension lag or on top of the finger to act as a PIP extension assist. It is most frequently used to prevent abduction or scissoring of fingers.*



SIRIS™ MEASURING INSTRUCTIONS

SILVER RING SPLINT COMPANY | PHONE: 434.971.4052 or 800.311.7028 | FAX: 434.971.8828 or 888.456.8828

CustomerService@SilverRingSplint.com | www.SilverRingSplint.com