



# SIRIS PULLEY RING SPLINT (P)

*Bowstringing of A-2 Pulley*


- 1** Choose a Round Ringsizer that fits at the base of the finger (MP crease).



- 2** Specify band width: 4, 6 or 10mm.

*Note: Wider band widths are tighter going over the joint. Splint will not be effective if the PIP joint is enlarged.*

*(Actual size shown)*



---

## SIRIS™ MEASURING INSTRUCTIONS

SILVER RING SPLINT COMPANY | PHONE: 434.971.4052 or 800.311.7028 | FAX: 434.971.8828 or 888.456.8828

*CustomerService@SilverRingSplint.com | www.SilverRingSplint.com*